



Anti- Flunitrazepam Polyclonal antibody (DPATB-H83102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Flunitrazepam Polyclonal Antibody
Target	Flunitrazepam
Immunogen	Synthetic peptide, 8amino derivative (part of the native molecule) conjugated to KLH via a diazo-coupling.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	250 µl
Buffer	Whole serum
Preservative	0.09% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Flunitrazepam, a benzodiazepine derivative, is a powerful hypnotic drug. It is also known as Rohypnol and sometimes used as a date rape drug. It has powerful hypnotic, sedative, anxiolytic, and skeletal muscle relaxant properties.

Keywords

Rohypnol;
