



Anti-Methandriol polyclonal antibody (DPATB-H83105)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Methandriol Polyclonal Antibody
Target	Methandriol
Immunogen	Methandriol (3beta)-BTG.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	100 µl
Buffer	20mM PBS, 150mM Sodium chloride
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Methandriol is the brand name for the anabolic steroid methylandrostenediol. Methandriol (or MAD) is one of the more rare and exotic anabolic steroids. It is actually 6androstenediol (5AD) that has had its chemical structure modified by adding a methyl group.

Keywords

MAD;