



Anti-Meperidine polyclonal antibody (DPATB-H83107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Meperidine Polyclonal Antibody
Target	Meperidine
Immunogen	Meperidine conjugated to Bovine thyroglobulin
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Bovine
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	1 mg, 10 mg
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Meperidine, (ethyl 1-methyl-4-phenylpiperidine-4-carboxylate), is a narcotic analgesic that can be used for the relief of most types of moderate to severe pain, including postoperative pain and the pain of labour. Prolonged use may lead to dependence of t

Keywords

Demerol; ethyl 1 methyl5phenylpiperidine5carboxylate; Isonipeccaine; Meperidol; Pethanol; Pethidin; Pethidine; Pethidimeter; Phetidine; Pipersal;