



Anti-GABA polyclonal antibody (DPATB-H83114)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-GABA Polyclonal Antibody
Target	GABA
Immunogen	Chemical / Small Molecule by a Glutaraldehyde linker.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Cellular Localization	Cell Membrane, Cytoplasmic, Plasma membrane and Secreted
Format	Liquid
Size	50 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Gamma-aminobutyric acid (GABA) is a major inhibitory neurotransmitter. GABA acts at inhibitory synapses in the brain and spinal cord. Inhibition is provoked by GABA binding resulting in hyperpolarization of the synaptic transmembrane potential of the affe

Keywords

4 aminobutanoic acid; GABA; Gamma amino butyric acid; Gamma aminobutyric acid;
