



Sheep anti-cGMP Polyclonal antibody (DPATB-H83142)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-cGMP Polyclonal Antibody
Target	cGMP
Immunogen	Succinylated cGMP conjugated to bovine serum albumin.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Bovine
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, RIA
Format	Liquid
Size	100 µl
Buffer	serum with 0.01% sodium azide
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Cyclic guanosine monophosphate (cGMP) serves as a second messenger in a manner similar to that observed with cAMP. Peptide hormones, such as the natriuretic factors, activate receptors that are associated with membrane-bound guanylate cyclase (GC). Recept

Keywords

Cyclic GMP; Cyclic guanosine monophosphate; Guanosine 36Cyclic Monophosphate;
