



# Anti-Choline Glutaric acid polyclonal antibody (DPATB-H83144)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Conjugated Choline Glutaric acid (conjugated) Polyclonal Antibody
Target	Conjugated Choline Glutaric acid
Immunogen	Synthetic choline glutaric anhydride (GA)-Carriers.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, ICC, IHC-P, FC
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C.

## BACKGROUND

**Introduction**

Choline increases carnitine uptake and retention. Choline supplements are inexpensive, safe (probably even in all children requiring anticholinergics) and can provide spectacular evidence of the suboptimal efficiency of carnitine supplementation by increa

---

**Keywords**

ACh;

---