



Anti-M. tuberculosis Polyclonal antibody (DPATB-H83149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Mycobacterium tuberculosis Polyclonal Antibody
Target	M. tuberculosis
Immunogen	Purified Protein Derivative
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	M. tuberculosis
Purification	IgG fraction
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr
Positive Control	Infected lung
Format	Liquid
Size	500 µl
Buffer	PBS, Carrier protein, Da Vinci Green Diluent, pH 7.3
Preservative	0.1% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Mycobacterium tuberculosis is the most common cause of tuberculosis. Primary infection begins with inhalation of 1 to 2 aerosolised bacilli. The pathogenicity of the organism is determined by its ability to escape host immune responses as well as elicit??

Keywords

M. tuberculosis;
