



# Anti-HSV type 1 Polyclonal antibody (DPATB-H83158)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-HSV1, prediluted Polyclonal Antibody
<b>Target</b>	HSV type 1
<b>Immunogen</b>	Whole rabbit corneal cells, infected with herpes simplex virus type 1 (strain Mac Intyre), were solubilized with detergent.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	HSV
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, IHC-Fr
<b>Format</b>	Liquid
<b>Buffer</b>	Inert stabilizer
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Store at +4°C. Do not freeze.

## BACKGROUND

**Introduction** Herpes simplex type 1 (HSV-1) belongs to a family that includes HSV-2, Epstein-Barr virus

(EBV) and Varicella zoster (chicken pox) virus amongst others. HSV-1 and HSV-3 are extremely difficult to distinguish from each other. Members of this family have a?

---

**Keywords**

Herpes simplex virus 1; Herpes Simplex Virus; Herpes Simplex Virus Type 1; HSV 1; Human herpesvirus 1; Human herpesvirus type 1;

---