



Anti-HSV type 2 Polyclonal antibody (DPATB-H83163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-HSV2, prediluted Polyclonal Antibody
Target	HSV type 2
Immunogen	Whole rabbit corneal cells, infected with herpes simplex virus type II (strain M5), were solubilized with detergent.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr
Format	Liquid
Buffer	Inert stabilizer
Preservative	0.05% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction Herpes simplex type3(HSV2) belongs to a family that includes HSV1, Epstein-Barr virus (EBV)

and Varicella zoster (chicken pox) virus. HSV1 and HSV3 are extremely difficult to distinguish from each other. These viruses have a DNA genome, an icosahedral p??

Keywords

Herpes Simplex Virus Type 2;
