



Anti-HSV type 1, 2 Polyclonal antibody (DPATB-H83164)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-HSV1 + HSV2, prediluted Polyclonal Antibody
Target	HSV type 1, 2
Immunogen	Whole rabbit corneal cells, infected with herpes simplex virus type I (strain Mac Intyre) and type II (strain M5), were solubilized with detergent.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-Fr, IHC-P
Format	Liquid
Buffer	Inert stabilizer
Preservative	0.05% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction There are two main types of Herpes Simplex Virus (HSV), 1 and 2. The primary difference

between the two types is in where they typically establish latency in the body. HSV1 usually establishes latency in the trigeminal ganglion, a collection of nerve cell

Keywords

Herpes Simplex Virus 1 and 2; Herpes simplex virus 1; Herpes simplex virus 2; Herpes simplex virus type 1; Herpes simplex virus type 2; HSV 1; HSV 2;
