



Anti-HSV type 1 Polyclonal antibody (DPATB-H83168)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-HSV1 Polyclonal Antibody
Target	HSV type 1
Immunogen	Whole rabbit corneal cells, infected with herpes simplex virus type I (strain Mac Intyre), were solubilized with detergent.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr, ICC, IF, ELISA
Format	Liquid
Size	500 µl
Buffer	1% BSA
Preservative	0.05% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Herpes simplex type 1 (HSV-1) belongs to a family that includes HSV-2, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus amongst others. HSV-1 and HSV-10 are extremely difficult to distinguish from each other. Members of this family have a

Keywords

Herpes simplex virus 1; Herpes Simplex Virus; Herpes Simplex Virus Type 1; HSV 1; Human herpesvirus 1; Human herpesvirus type 1;
