



Anti-T. gondii Polyclonal antibody (DPATB-H83199)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Toxoplasma gondii Polyclonal Antibody
Target	T. gondii
Immunogen	Toxoplasma gondii
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. gondii
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-P
Positive Control	Infected brain.
Format	Liquid
Size	500 µl
Buffer	1% BSA, 10mM PBS, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Toxoplasma is a crescent shaped sporozoan that lives as an intracellular parasite in various tissues of many vertebrates and completes its life cycle in a single host. Its life cycle includes two phases called the intestinal (or enteroepithelial) and extr

Keywords

T. gondii,
