



## Anti-T. gondii SAG1 protein monoclonal antibody (DPATB-H83211)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	SAG1 Protein
<b>Immunogen</b>	T. gondii extract
<b>Isotype</b>	IgG2a
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	T. gondii
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(cap), WB, IF
<b>Format</b>	Purified, Liquid
<b>Size</b>	500 µg
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.