



Anti-S. pneumoniae type 18 Polyclonal antibody (DPATB-H83212)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Guinea Anti-Streptococcus pneumoniae type29Polyclonal Antibody
Target	S. pneumoniae type 18
Immunogen	Extract of species specific membrane proteins.
Isotype	IgG
Source/Host	Guinea pig
Species Reactivity	S. pneumoniae
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at-20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Bronchial pneumonia is most prevalent in infants, young children and aged adults. It is caused by various bacteria, including <i>Streptococcus pneumoniae</i> . <i>Streptococcus pneumoniae</i> are Gram positive, lancet shaped cocci (elongated cocci with a slightly pointed
Keywords	Pneumococcal surface protein A,pspA,