



Anti-Staphylococcus Group B Polyclonal antibody (DPATB-H83213)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus Group B Polyclonal Antibody
Target	Staphylococcus Group B
Immunogen	Streptococci, group B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Staphylococcus
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ICC/IF
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2. No stabilizing proteins have been added.
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Streptococci form part of the normal human flora that resides on the skin, and can also colonise the respiratory, gastrointestinal, and genitourinary tracts. Streptococci can cause a range of diseases, from the less serious but common sore throats and ski

Keywords Beta strep, Group b streptococci, Group b streptococcus, Strep b,
