



# Anti-T. pallidum Polyclonal antibody (DPATB-H83214)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Treponema pallidum Polyclonal Antibody
Target	T. pallidum
Immunogen	Highly purified Treponema pallidum.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. pallidum
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Treponema pallidum is the causative agent of syphilis. It is a gram negative spirochete, a helical to sinusoidal bacterium with outer and cytoplasmic membranes, a thin peptidoglycan layer, and periplasmic flagella.

---

**Keywords**

T. pallidum,

---