



Anti-T. gondii Polyclonal antibody (DPATB-H83215)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Toxoplasma gondii Polyclonal Antibody
Target	T. gondii
Immunogen	RH strain of Toxoplasma gondii purified Trophozoites.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. gondii
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, ICC/IF
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Toxoplasma is a crescent shaped sporozoan that lives as an intracellular parasite in various tissues of many vertebrates and completes its life cycle in a single host. Its life cycle includes two phases called the intestinal (or enteroepithelial) and extr

Keywords

T. gondii,
