



Anti-Serotonin polyclonal antibody (DPATB-H83240)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Serotonin Polyclonal Antibody
Target	Serotonin
Immunogen	Serotonin conjugated to bovine thyroglobulin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-Fr, IHC-P
Cellular Localization	Cytoplasmic and Secreted
Positive Control	Small Intestine
Format	Liquid
Size	500 µl
Buffer	Whole serum diluted with PBS, 1% BSA.
Preservative	0.09% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Serotonin (5-hydroxytryptamine, or 5-HT) is a monoamine neurotransmitter synthesised in serotonergic neurons in the central nervous system and enterochromaffin cells in the gastrointestinal tract. Serotonin plays an important part in the biochemistry of d
Keywords	5 HT, 5 hydroxytryptamine, 5HT,