



Anti-Quinaldic Acid polyclonal antibody (DPATB-H83253)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Quinaldic Acid Polyclonal Antibody
Target	Quinaldic Acid
Immunogen	Synthetic peptide Quinaldic Acid conjugated to a protein carrier.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The N-heterocycle quinaldic acid (quinoline 2-carboxylic acid) is utilized by *Alcaligenes* sp. F-2

as sole source of carbon and energy.

Keywords

Quinoline3carboxylic acid,
