



Anti-Vitamin B3 polyclonal antibody (DPATB-H83259)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Vitamin B3 Polyclonal Antibody
Target	Vitamin B3
Immunogen	Chemical / Small Molecule, Vitamin B3 conjugated to BSA.
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Preadsorbed IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Vitamin B3 is a water-soluble vitamin whose derivatives such as NADH, NAD, NAD⁺, and NADP play essential roles in energy metabolism in the living cell and DNA repair. The designation Vitamin B3 also includes the amide form, nicotinamide or niacinamide. Se

Keywords

Niacin, Nicotinamide, Nicotinic acid,
