



Anti-Thioctic Acid polyclonal antibody (DPATB-H83261)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Thioctic acid Polyclonal Antibody
Target	Thioctic acid
Immunogen	Thioctic acid conjugated to BSA.
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Preadsorbed IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Thioctic acid is a natural anti-oxidant which has been used for years in Europe because of its ability to help maintain or restore liver health. Its strong anti-oxidant capabilities have been used for radiation sickness and diabetes-related neuropathy.

Keywords

Alpha lipoic acid,
