



Anti-P. mirabilis Polyclonal antibody (DPATB-H83266)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Proteus mirabilis Polyclonal Antibody
Target	P. mirabilis
Immunogen	Proteus mirabilis total proteins.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	P. mirabilis
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Proteus mirabilis is a Gram-negative, facultatively anaerobic bacterium. It shows swarming, motility, and urease activity. P. mirabilis causes 90% of all 'Proteus' infections. It belongs to the Tribe Proteae.

Keywords

p. mirabilis, Proteus mirabilis,
