



Anti-Tyrosine polyclonal antibody (DPATB-H83268)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Tyrosine Polyclonal Antibody
Target	Tyrosine
Immunogen	Synthetic tyrosine conjugated to BSA by a Glutaraldehyde linker.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Preadsorbed antiserum
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Tyrosine (4-hydroxyphenylalanine, or 2-amino-3(4-hydroxyphenyl)-propanoic acid) is one of the 22 amino acids that are used by cells to synthesize proteins. Tyrosine cannot be completely synthesized by animals, although it can be made by hydroxylation of p

Keywords

2 amino 3(4 hydroxyphenyl) propanoic acid, 4 hydroxyphenylalanine,