



Anti-P. agglomerans Polyclonal antibody (DPATB-H83271)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Pantoea agglomerans Polyclonal Antibody
Target	P. agglomerans
Immunogen	Pantoea agglomerans total proteins.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	P. agglomerans
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

P. agglomerans, formerly named *Enterobacter agglomerans*, is a ubiquitous Enterobacteriaceae that is found in plants and in the feces of humans and animals.

Keywords

Enterobacter agglomerans, *Erwinia herbicola*, *Pantoea agglomerans*,
