



Anti-Rotenone polyclonal antibody (DPATB-H83272)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Rotenone Polyclonal Antibody
Target	Rotenone
Immunogen	Chemical / Small Molecule: Rotenone conjugated to a protein carrier.
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Rotenone is a colorless-to-red, odorless solid. In solution it is used as a broad spectrum
---------------------	--

insecticide that works by inhibiting the transfer of electrons from Fe-S centers in Complex I to ubiquinone (see electron transfer chain). This prevents NADH from b

Keywords

Canex,
