



Anti-Vitamin B1 polyclonal antibody (DPATB-H83275)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Vitamin B1 Polyclonal Antibody
Target	Vitamin B1
Immunogen	Synthetic vitamin B1 conjugated to BSA
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Preadsorbed antiserum
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Also known as thiamine, vitamin B1 plays an important role in helping the body convert carbohydrates and fat into energy. It is essential for normal growth and development and helps to maintain proper functioning of the heart and the nervous and digestive

Keywords

Aberic acid,Thiamin,Thiamine,Vitamin B1,
