



# Anti-Tetanus Toxin polyclonal antibody (DPATB-H83287)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Tetanus Toxin Polyclonal Antibody
Target	Tetanus Toxin
Immunogen	Tetanus toxoid (formaldehyde inactivated tetanus toxin).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Clostridium tetani
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C.

## BACKGROUND

<b>Introduction</b>	Tetanus toxin acts by inhibiting neurotransmitter release. It binds to peripheral neuronal synapses, is internalized and moves by retrograde transport up the axon into the spinal cord where it can move between postsynaptic and presynaptic neurons. It inhi
<b>Keywords</b>	Tentoxylisin,Tetanus toxin chain H,Tetanus toxin chain L,Tetanus toxin heavy chain,Tetanus toxin light chain,