



Anti-Thiabendazole polyclonal antibody (DPATB-H83289)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Thiabendazole Polyclonal Antibody
Target	Thiabendazole
Immunogen	Thiabendazole conjugated to Bovine thyroglobulin
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Bovine
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Thiabendazole is a 2-Substituted benzimidazole, 2-(1,3-thiazol-4-yl)-1H-benzimidazole, and is active against a variety of nematodes. It is specifically used for strongyloidiasis. Thiabendazole has CNS side effects and hepatotoxic potential.

Keywords

2 (1,3 thiazol5yl) 1H
benzimidazole,Bioguard,Bovizole,Equizole,Lombristop,Mintezol,Thiabendazol,Thiabenzole,
