



Anti-Ritodrine polyclonal antibody (DPATB-H83290)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Ritodrine Polyclonal Antibody
Target	Ritodrine
Immunogen	Ritodrine conjugated to BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Ritodrine, (4-[(1R,2S)-1-hydroxy-2-[2-(4-hydroxyphenyl)ethylamino]propyl]phenol), is an adrenergic beta-agonist used to control premature labour. It has a molecular weight of 287.35354 g/mol and molecular formula C₁₇H₂₁NO₃.

Keywords

4 [(1R,2S) 1 hydroxy3[2 (4 hydroxyphenyl)ethylamino]propyl]phenol, Ritodrina, Ritodrinium, Yutopar,