



Anti-S. maltophilia Polyclonal antibody (DPATB-H83322)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Stenotrophomonas maltophilia Polyclonal Antibody
Target	S. maltophilia
Immunogen	Preparation of Stenotrophomonas maltophilia total proteins
Isotype	IgG
Source/Host	Rat
Species Reactivity	S. maltophilia
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Stenotrophomonas maltophilia (previously known as Pseudomonas maltophilia) is a Gram negative bacterium that is commonly found in a variety of aquatic environments. It is intrinsically multi-drug resistant and occasionally causes bacteraemic and organ-spe

Keywords

Pseudomonas maltophilia, S. maltophilia, SM, Stenotrophomonas maltophilia,