



Anti-L. monocytogenes Polyclonal antibody (DPATB-H83323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Listeria monocytogenes Polyclonal Antibody
Target	L. monocytogenes
Immunogen	Cell preparation: Listeria monocytogenes
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	L. monocytogenes
Purification	Protein A purified
Conjugate	Unconjugated
Applications	Conjugation, ELISA
Format	Liquid
Size	100 µl
Buffer	10mM PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

The genus *Listeria* comprises six species: *L. monocytogenes*, *L. innocua*, *L. welshimeri*, *L. seeligeri*, *L. ivanovii* and *L. grayi*. *Listeria monocytogenes*, the most commonly isolated pathogenic member, is associated with a wide spectrum of human and animal dis

Keywords

L. monocytogenes;
