



# Anti-E. coli H-28 Polyclonal antibody (DPATB-H83329)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-E. coli H-28 Polyclonal Antibody
<b>Target</b>	E. coli H-28
<b>Immunogen</b>	Tissue / cell preparation of E. coli H-28 (heat inactivated).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	E. coli
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Agg
<b>Positive Control</b>	Antigenic suspension.
<b>Format</b>	Liquid
<b>Size</b>	5 ml
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C. Do not freeze.

## BACKGROUND

**Introduction**

Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae - bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es

---

**Keywords**

Escherichia coli H-28;

---