



# Anti-E. coli O151 Polyclonal antibody (DPATB-H83340)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli O151 Polyclonal Antibody
Target	E. coli O151
Immunogen	Tissue / cell preparation of E. coli O151 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

## BACKGROUND

<b>Introduction</b>	Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type)
<b>Keywords</b>	Escherichia coli O151;