



# Anti-E. coli O153 Polyclonal antibody (DPATB-H83341)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                               |
|---------------------------|---------------------------------------------------------------|
| <b>Product Overview</b>   | Rabbit Anti-E. coli O153 Polyclonal Antibody                  |
| <b>Target</b>             | E. coli O153                                                  |
| <b>Immunogen</b>          | Tissue / cell preparation of E. coli O153 (heat inactivated). |
| <b>Isotype</b>            | IgG                                                           |
| <b>Source/Host</b>        | Rabbit                                                        |
| <b>Species Reactivity</b> | E. coli                                                       |
| <b>Purification</b>       | Whole antiserum                                               |
| <b>Conjugate</b>          | Unconjugated                                                  |
| <b>Applications</b>       | Agg                                                           |
| <b>Positive Control</b>   | Antigenic suspension.                                         |
| <b>Format</b>             | Liquid                                                        |
| <b>Size</b>               | 2 ml                                                          |
| <b>Buffer</b>             | Whole serum                                                   |
| <b>Preservative</b>       | None                                                          |
| <b>Storage</b>            | Store at +4°C. Do not freeze.                                 |

## BACKGROUND

**Introduction** Escherichia coli is a group of gram-negative bacilli belonging to the family Enterobacteriaceae and is considered as one of the bacteria in the human intestinal normal microbial flora. Their serological types are determined in combination with somatic antigens.

---

**Keywords** Escherichia coli O153;

---