



Anti-E. coli O157 Polyclonal antibody (DPATB-H83342)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli O157 Polyclonal Antibody
Target	E. coli O157
Immunogen	Tissue / cell preparation of E. coli O157 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type)

Keywords Escherichia coli O157;
