



# Anti-E. coli O28ac Polyclonal antibody (DPATB-H83348)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                |
|---------------------------|----------------------------------------------------------------|
| <b>Product Overview</b>   | Rabbit Anti-E. coli O28ac Polyclonal Antibody                  |
| <b>Target</b>             | E. coli O28ac                                                  |
| <b>Immunogen</b>          | Tissue / cell preparation of E. coli O28ac (heat inactivated). |
| <b>Isotype</b>            | IgG                                                            |
| <b>Source/Host</b>        | Rabbit                                                         |
| <b>Species Reactivity</b> | E. coli                                                        |
| <b>Purification</b>       | Whole antiserum                                                |
| <b>Conjugate</b>          | Unconjugated                                                   |
| <b>Applications</b>       | Agg                                                            |
| <b>Positive Control</b>   | Antigenic suspension.                                          |
| <b>Format</b>             | Liquid                                                         |
| <b>Size</b>               | 2 ml                                                           |
| <b>Buffer</b>             | Whole serum                                                    |
| <b>Preservative</b>       | None                                                           |
| <b>Storage</b>            | Store at +4°C. Do not freeze.                                  |

## BACKGROUND

**Introduction** Escherichia coli is a group of gram-negative bacilli belonging to the family Enterobacteriaceae and is considered as one of the bacteria in the human intestinal normal microbial flora. Their serological types are determined in combination with somatic antigens.

---

**Keywords** Escherichia coli O28ac;

---