



Anti-E. coli pili 987P Polyclonal antibody (DPATB-H83355)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli pili 987P Polyclonal Antibody
Target	E. coli pili 987P
Immunogen	Tissue / cell preparation of toxigenic E.coli pili 987P (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	5 ml
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae, bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es

Keywords

Escherichia coli pili 987P;