



Anti-E. coli Polyclonal antibody (DPATB-H83359)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli Polyclonal Antibody
Target	E. coli
Immunogen	Tissue / cell preparation (Mixture of E. coli serotypes: K12, O111, O55, O125, O20, and O157).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Size	500 µl
Buffer	pH: 7.20Preservative: 0.1% Sodium azideConstituent: PBS
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es

Keywords

Escherichia coli;
