



Anti-E. coli O157 Polyclonal antibody (DPATB-H83366)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-E. coli O157 Polyclonal Antibody
Target	E. coli O157
Immunogen	Tissue / cell preparation (whole cells of Escherichia coli O157:H7)
Isotype	IgG
Source/Host	Goat
Species Reactivity	E. coli
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µg
Buffer	40% Glycerol, PBS
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O175) and flagella antigens (H type)

Keywords Escherichia coli;
