



Anti-E. coli Polyclonal antibody (DPATB-H83374)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-E. coli-Azide free Polyclonal Antibody
Target	E. coli
Immunogen	Heat interacted blend of 7 different strains of Escherichia coli : TOP 10F, HB101, JM109, BL21, DH5 alpha and K12.
Isotype	IgG
Source/Host	Goat
Species Reactivity	E. coli
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	500 µl
Buffer	PBS, Sodium citrate, pH 7.0-7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es

Keywords

Escherichia coli;
