



Anti-E. coli K88AB pilus Polyclonal antibody (DPATB-H83379)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-E. coli K88AB pilus Polyclonal Antibody
Target	E. coli K88AB pilus
Immunogen	E. coli strain 08:K87:K88ab
Isotype	IgG
Source/Host	Sheep
Species Reactivity	E. coli
Purification	Ammonium Sulphate Precipitation
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Cell Membrane
Format	Liquid
Size	100 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Escherichia coli is the head of the large bacterial family, Enterobacteriaceae, the enteric bacteria, which are facultatively anaerobic Gram negative rods that live in the intestinal tracts of animals in health and disease. Pili are macromolecular structur
Keywords	Escherichia coli; Escherichia coli K88AB pilus;