



Anti-Neurofibrillary Tangles polyclonal antibody (DPATB-H83394)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Neurofibrillary tangles Polyclonal Antibody
Target	Neurofibrillary tangles
Immunogen	Disaggregated paired helical filaments.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P, ELISA, WB
Positive Control	Alzheimer disease brain.
Format	Liquid
Size	50 µl
Buffer	Whole Serum
Preservative	0.09% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Neurofibrillary tangles consist of insoluble twisted fibers that are found inside of the brain's cells. They primarily consist of a protein called tau, which forms part of a structure called a microtubule. The microtubule helps transport nutrients and oth
Keywords	NFT;