



Anti-Neurofibrillary Tangles polyclonal antibody (DPATB-H83396)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Neurofibrillary tangles Polyclonal Antibody
Target	Neurofibrillary tangles
Immunogen	Neurofibrillary tangles which were extracted in boiling SDS and purified by sucrose gradient centrifugation.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr, ELISA
Positive Control	Human Alzheimers brain tissue
Format	Liquid
Size	100 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Neurofibrillary tangles consist of insoluble twisted fibers that are found inside of the brain's cells. They primarily consist of a protein called tau, which forms part of a structure called a microtubule. The microtubule helps transport nutrients and oth

Keywords

NFT;
