



Anti-Diethylstilbestrol polyclonal antibody (DPATB-H83408)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Diethylstilbestrol Polyclonal Antibody
Target	Diethylstilbestrol
Immunogen	Diethylstilbestrol conjugated to KLH by a hemisuccinate linker.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Diethylstilbestrol (DES) is a synthetic form of the hormone estrogen that was prescribed to pregnant women between about 1940 and 1971 because it was thought to prevent miscarriages. DES may increase the risk of uterine, ovarian, or breast cancer in women

Keywords

DES;
