



Anti-serum polyclonal antibody (DPATB-H83410)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Human serum (C4 free) Polyclonal Antibody
Target	Human serum
Immunogen	Human serum, C4 free
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, WB, IEP, ID
Format	Liquid
Size	1 ml
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

The human serum proteome as a whole reflects changes in the physiological state of the human body. The onset of diseases invariably induces changes in the protein content, perhaps by producing more of one protein, less of another, or by producing a whole

Keywords

C4 Anaphylatoxin; C4A2; C4A3; C4A4; C4A6; C4B; C4S; CO4; Rodgers Form Of C4;
