



Magic™ Anti-Hirudin (full length) polyclonal antibody (DPATB-H83415)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Hirudin Polyclonal Antibody
Target	Hirudin
Immunogen	Purified full length recombinant Hirudin (Hirudo medicinalis), a 65 aa polypeptide.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	IgY fraction
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µl
Buffer	0.9% Sodium chloride
Preservative	0.1% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Hirudin is a highly effective and very specific anticoagulant which specifically inhibits formation of fibrin and activation of factors V, VIII and XIII by forming irreversible, noncovalent 1:1 complexes with a-thrombin. Hirudin is a 65 amino acid polypep

Keywords

Hirudin 1; Lepirudin; Repludin;
