



Anti-Thymulin polyclonal antibody (DPATB-H83416)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Thymulin Polyclonal Antibody
Target	Thymulin
Immunogen	Synthetic Thymulin EAKSQGGSN (Human). EAKSQGGSN Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, RIA
Cellular Localization	Secreted
Format	Liquid
Size	50 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Thymulin, a nonapeptide secreted by the thymus, is an essential hormone for T lymphocyte differentiation and function.

Keywords

Serum thymic factor,
