



# Anti-T. gondii Polyclonal antibody (DPATB-H83417)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Toxoplasma gondii, prediluted Polyclonal Antibody
Target	T. gondii
Immunogen	Toxoplasma gondii
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. gondii
Conjugate	Unconjugated
Applications	IHC-P
Positive Control	Infected brain tissue.
Format	Prediluted
Preservative	None
Storage	Store at +4°C.

## BACKGROUND

Introduction	Toxoplasma is a crescent shaped sporozoan that lives as an intracellular parasite in various tissues of many vertebrates and completes its life cycle in a single host. Its life cycle includes two phases called the intestinal (or enteroepithelial) and extr
--------------	---

**Keywords**

T. gondii,

---