



Anti-Triiodothyronine polyclonal antibody (DPATB-H83419)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Triiodothyronine-Azide free Polyclonal Antibody
Target	Triiodothyronine
Immunogen	T3-BgG
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A or G purified
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µl
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Preservative	0.05% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	A hormone made by the thyroid gland. It has three iodine molecules attached to its molecular structure. It is the most powerful thyroid hormone, and affects almost every process in the body, including body temperature, growth, and heart rate.
Keywords	liothyronine, T3,