



Anti-Triiodothyronine polyclonal antibody (DPATB-H83420)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Triiodothyronine Polyclonal Antibody
Target	Triiodothyronine
Immunogen	T3 conjugated to BSA
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	A hormone made by the thyroid gland. It has three iodine molecules attached to its molecular structure. It is the most powerful thyroid hormone, and affects almost every process in the body, including body temperature, growth, and heart rate.
Keywords	liothyronine,T3,